


# Increase Testing Efficiency

	 <b>Q-global</b>	 <b>Q-interactive</b>
<b>Test Administration</b>	<p>The following are available for select assessments</p> <ul style="list-style-type: none"> <li>On-screen assessments (client enters responses directly into digital device)</li> <li>Remote on-screen administration (a link is sent to respondent for completion of assessment)</li> <li>Group administration</li> <li>Paper-and-pencil administration with raw scores input into Q-global for scoring.</li> </ul>	<p>No more lugging heavy test kits, just two iPads.</p> <ul style="list-style-type: none"> <li>Batteries can be selected, tailored, or adjusted on-the-fly</li> <li>The clinician controls the administration of the test</li> <li>Test instructions, stimuli, record forms, norms tables are all embedded in the system and accessible from the clinician's iPad</li> <li>Responses are captured via handwriting and audio recording.</li> </ul>
<b>Scoring</b>	<ul style="list-style-type: none"> <li>Input raw scores from paper and pencil</li> <li>Administrations into a computer at the conclusion of the battery</li> <li>At the end of the test and with Internet connectivity, view test results.</li> </ul>	<ul style="list-style-type: none"> <li>Monitor test performance with on-the-fly scoring</li> <li>Have the flexibility to override scores when the client self corrects</li> <li>No need for internet connectivity to view score results.</li> </ul>
<b>Reporting</b>	<ul style="list-style-type: none"> <li>Basic score reports, and combination reports</li> <li>interpretive reports for selected products are available.</li> </ul>	<ul style="list-style-type: none"> <li>Online scoring reports are available</li> <li>Digital notes, which are captured directly on the iPad, can be incorporated into reports.</li> </ul>

For more information about Q-global visit:  
[pearsonclinical.com.au/QGlobal](http://pearsonclinical.com.au/QGlobal)



**recycle**  
 When you have finished with this batch please recycle it.  
 Pearson Education Ltd is committed to reducing its impact on the environment by using responsibly sourced recycled paper.

# Q-global® and Q-interactive® – providing flexible digital options

You can choose either or both options, to maximise your digital assessment experience.

## System-Driven Take On-Screen Assessment



- Supervised, group or remote administration
- The clinician is not required during testing session
- The client responds directly on the digital device
- Test administration can be completed anywhere internet is available
- Pay per report

## Generate Scores & Reports from Paper Test


- Clinician administers test in paper-and-pencil and captures raw scores on the record form
- Raw scores are entered into the computer
- Score and interpretive reports are automatically generated
- Pay per report or unlimited scoring subscription.



## Clinician-Driven Administer and Score Individual Assessments Using Two iPads



- Clinician administers tests directly to the client
- Scoring is on-the-fly
- Central** is the secured, web-based application and where a clinician manages client profiles, creates test batteries, runs reports, and exports/archives scored batteries.
- Assess** is the secured, tablet-based Q-interactive application where a clinician can administer a test using two iPads even without an internet connection.



For more information about Q-interactive visit:  
[pearsonclinical.com.au/qi](http://pearsonclinical.com.au/qi)

